

CoreFlow DOWNLOAD AND INSTALLATION INSTRUCTIONS

The latest stable release of CoreFlow is available for download at
<http://pawsonlab.mshri.on.ca/CoreFlow/download.php>.

Installation Instructions

INSTALLATION REQUIREMENTS.....	1
INSTALLING CoreFlow FOR THE FIRST TIME.....	1
UPGRADING TO A NEWER VERSION OF CoreFlow.....	3

INSTALLATION REQUIREMENTS

- ✓ Apache HTTP server: <http://httpd.apache.org/> - on both web and application servers
- ✓ PHP version 5 or later: <http://php.net/> with MySQL support (php-mysql) - web and application servers
- ✓ MySQL DBMS version 5 or later: <http://www.mysql.com/> - web and application servers
- ✓ Perl: <http://perl.org> with the following modules: - web and application servers
 - DBI
 - DBD::mysql
 - XML::XPath
 - XML::Parser
 - If required by above XML modules, install Expat XML parser: <http://expat.sourceforge.net/>
- ✓ R statistical software: <http://www.r-project.org/> - application server
- ✓ Python version 2.4 or later: <http://python.org/> with MySQLdb module installed - application server
- ✓ BioPython: http://biopython.org/wiki/Main_Page - application server
- ✓ Optional: Install git from <http://git-scm.com/>. For more information on the use of git for management of your analysis scripts developed in task playgrounds please refer to the YouTube videos on <http://pawsonlab.mshri.on.ca/CoreFlow>.
- ✓ **The application server must be configured to allow execution of shell scripts from the runR.pl

INSTALLATION PROCEDURE

INSTALLING CoreFlow FOR THE FIRST TIME

- CoreFlow can be installed on any Unix/Linux/Mac platform
 - Database and server code can be installed either on the same hardware node, or on separate nodes, according to your preference.
1. Download and extract [coreFlow_web_2013_01.tgz](#) into your web server directory.

The unpacked directory contains:

- A folder named CoreFlow – main website directory

- A database dump file, named **mol_bio_dump.sql**, and
 - A copy of the Installation Instructions in PDF format
2. Copy the R subfolder from the unpacked CoreFlow directory onto your application server. The application server will run the R, Perl, Python/BioPython, etc. scripts that are attached in each data analysis task.
 3. Define the storage location of R CGI files on the (possibly different) application server:

In **.config.xml**, edit the following section:

```
<R>
  <cgi_dir>http://your\_app\_server\_name/R\_directory/</cgi_dir>      <!-- e.g. This may be
  "http://localhost/CoreFlow/R/" -->
  <temp_dir>full_path_to_any_temp_directory</temp_dir>      <!-- e.g.
  /tmp/CoreFlow/tmp; if this directory does not exist, create it, AND change permissions to
  a+rwX -->
</R>
```

4. Edit your Apache httpd.conf file as follows (if the web server and application server are on different machines, perform the configuration on both machines):

- Define ServerName, DocumentRoot and ServerAlias for your CoreFlow setup directory.
- Set LD_LIBRARY_PATH to point to your MySQL directory.
- It is helpful to designate an ErrorLog file, e.g. /my_install_dir/**my_coreflow-error.log**.

Example:

(replace the highlighted directives according to your server's specifications)

```
<Directory /var/www/html/CoreFlow>
  IndexIgnore .config.xml
  DirectoryIndex index.php
  Options Indexes FollowSymLinks
  AddType application/x-httpd-php .php      # MUST HAVE THIS TO VIEW PHP PAGES!!!!
  Options ExecCGI
  AddHandler cgi-script .cgi .pl .py      # DO NOT INCLUDE .PHP HERE!!!
  AllowOverride AuthConfig
  Order allow,deny
  Allow from all
</Directory>
```

5. On your database server, create an empty MySQL schema that will be populated in the next step:

```
mysql> create database mol_bio;
```

(replace *mol_bio* with a schema name of your choice)

6. Use the database dump file **mol_bio.sql** to add tables to the newly created schema. At the command prompt, type:

```
xterm> mysql -u root -p mol_bio < mol_bio.sql
```

7. Create two new MySQL users – a web user, e.g. 'coreflow_web', and an app user, e.g. 'coreflow_app':

```
mysql> create user coreflow_web identified by 'web_passwd';
```

```
mysql> create user coreflow_app identified by 'app_passwd';
```

8. Grant SELECT, INSERT and UPDATE privileges on database *mol_bio* to the new users:

```
mysql> GRANT SELECT, CREATE TEMPORARY TABLES, LOCK TABLES ON mol_bio.* TO  
'coreflow_web'@'your_web_server_name_or_ip_address';
```

```
mysql> GRANT SELECT, INSERT, UPDATE, DELETE, ALTER, DROP, LOCK TABLES, CREATE  
TEMPORARY TABLES ON mol_bio.* TO 'coreflow_app'@'your_web_server_name_or_ip_address';
```

```
mysql> flush privileges;
```

9. In the CoreFlow folder, edit *.config.xml*:

```
<config>  
  <db> <!-- related to the database mysql server -->  
  <host>YOUR DATABASE SERVER HOSTNAME</host>  
  <database>mol_bio</database>  
  <user passwd="app_passwd">coreflow_app</user>  
  <user_public passwd="web_passwd">coreflow_web</user_public>  
  ...  
</config>
```

For users who are allowed to insert/delete/update (Write) into DB or eXecute files, add their workstations' IP addresses to the following line:

```
  <allow_wx_IPs>(user1\workstation\IP\address)|(user2\workstation\IP\address)|...  
</allow_wx_IPs>
```

THERE MUST BE ****NO SPACES**** IN THIS LINE!!!

10. Copy *.config.xml* from your CoreFlow web server directory into the R subdirectory on the application server

UPGRADING TO A NEWER VERSION OF CoreFlow

- 1) Make a copy of your current *.config.xml* file
- 2) Download the *.tgz* and the *mysql_dump* files from the 'Download' section of the public CoreFlow website
- 3) Erase your current CoreFlow installation directory and extract the downloaded archive
- 4) Replace the database tables with the latest downloaded mysql dump.
- 5) Copy back your saved *.config.xml* file into the extracted CoreFlow directory